

GTR: DSR Evaluation and Geometry Recommendation

**Asphalt Binder Expert Task Group
Salt Lake City, UT
April 27, 2016**

**Darin Hunter
Anton Paar USA
Ashland, VA**

Summary of Situation

Current Status

- ▶ **GTR modification of asphalt binder** provides environmental and performance benefits
- ▶ Current Superpave test methods (PG) are not suitable for rubber modified materials with **particle sizes up to 2+ mm** (mesh 10)
- ▶ Acceptance of new materials by highway agencies depends on suitable test methods for **rheological characterisation** (DSR)



Status

Scope of Work

- ▶ Evaluation of **various DSR measuring systems**
- ▶ Investigation of the **variability** in measured properties



Scope

Main Goals

- ▶ **Identify** practical and rheological **challenges**
- ▶ **Apply** relevant **engineering principles** such as measuring system standards
- ▶ **Suggest suitable measuring geometry** for PG measurements of GTR modified Asphalt Binder which provides the truest rheological data




Goals

Summary of Problem and Solution Approaches


Problem

- ▶ Swollen 2 mm rubber particles do not fit in the 1 or 2 mm gap if AASHTO T315
- ▶ Particle interferences impact results
- ▶ Accepted Governing Rules for Good Rheology
 - Measuring gap should be at least 5X and preferably 10X the largest particle size
 - AASHTO T315 / ASTM D7175 (DSR): max. 400 μm particle size (2mm gap, 8mm \varnothing plate)
- ▶ Questionable rheological data

1
mm



2
mm



Solution Approaches

- ▶ Increase the measuring gap
 - Plate-plate geometry with large gaps (up to 4mm)
 - Concentric cylinder geometry with large gaps
- ▶ Considerations:

<ul style="list-style-type: none"> • Does gap matter? • Does the data agree with conventional measuring gaps? 	<ul style="list-style-type: none"> • Are new sources of error introduced? • Are new practical concerns introduced?
---	--

4
mm






Plate-Plate Geometry: Pros and Cons

Pros

- 👍 Less sample volume
- 👍 Faster heating/cooling rates and temperature equilibrium
- 👍 Accustomed and approved system

Cons

- 👎 Reliable trimming and gap filling
 - Larger influence with smaller plates
- 👎 Representative sample volume
- 👎 Ratio diameter / gap → Geometry factors, Calibration (ISO 6712-10, DIN 53019-1)
- 👎 Time effects due to sagging (flow) of binder out of the measuring gap




Concerns

- ? Repeatability
- ? Reproducibility
- ? Reliability

**Sample: 10 % 10 mesh
GTR in PG64-22**

Plate-Plate Geometry: Sample Sag (Flow) Out of Gap


PP25/PE with 4 mm gap size - leak of binder over time



	$ G^* $ [Pa] $t = 0$ min	$ G^* $ [Pa] $t = 30$ min	$ G^* $ [Pa] $t = 60$ min	$ G^* $ [Pa] $t = 180$ min
Deviation referred to $t = 0$ min [%]	---	- 6.1	- 10.1	- 14.1


→ Flow of sample out of gap produces increasingly lower measured results over time even within time scale of 10 min thermal equilibrium

→ Parallel plate system should be used with maximum gap size of 2 mm




Remaining Rubber Matrix after 180 min Measurement (binder pooled on lower plate)

Plate-Plate Geometry: Repeatability



- ▶ Single, expert level user
- ▶ Issues contributing to lack of repeatability
 - Issues with; homogeneity, representative sample and trimming
 - Sample pulled out of gap and/or rubber particles at edge cannot be trimmed thus hang outside gap introducing additional drag
- ▶ New source of error introduced due to sampling and trimming issues

Plate-Plate Geometry: Reliable (Real) Measured Data



	Mean G^* at 46 °C [Pa]	Mean G^* at 16 °C [Pa]
PP08, 2 mm gap size	145537	19574667
PP08, 4 mm gap size	194573	24709667
Deviation with 4 mm gap size [%]	+ 33.69	+ 26.23

- ▶ Compared with data from the accepted 2mm gap
- ▶ No longer an absolute measuring geometry at large gap sizes
- ▶ Rheological data should be gap independent. When it ceases to be gap independent, the gap is too large.

Concentric Cylinder Geometry: Pros and Cons

Pros


- 👍 Large gaps possible
- 👍 Sample cannot flow out of the shear gap
 - No binder sag errors
- 👍 No trimming required
- 👍 Representative sample volume
- 👍 No edge effects
- 👍 Disposable cups can be used




Cons

- 👎 Increased sample volume
 - Slower heating and cooling rates
 - Longer temperature equilibrium time (at least 30 min)
- 👎 Relative measuring geometry requires calibration with traceable standard (Cannon N2700000SP)


Concentric Cylinder Geometry: Overview



Small Bob in Large Cup



Absolute (req. small gap)



Relative (large gap)

CC17SP


- Gap ≈ 6.2 mm
- Suitable for neat and RTFO asphalt binder (high PG temperatures)

CC10SP

- Gap ≈ 9.5 mm
- Suitable for PAV asphalt binder (intermediate PG temperatures)

Concentric Cylinder Geometries with Large Gap: Engineering Details

Concentric cylinder (CC) geometries with large gap (ratio $\delta_{cc} \gg 1.0847$)



B-CC17SP:

$$\delta_{cc} = \frac{R_e}{R_i} = \frac{14.46 \text{ mm}}{8.33 \text{ mm}} = 1.74 \quad (\text{Gap: } \approx 6.2 \text{ mm})$$

B-CC10SP:

$$\delta_{cc} = \frac{R_e}{R_i} = \frac{14.46 \text{ mm}}{5.00 \text{ mm}} = 2.89 \quad (\text{Gap: } \approx 9.5 \text{ mm})$$

- ▶ The combination of smaller cylinders (bobs) in a CC27 cup results in relative concentric cylinder systems (non ISO 3219 compliant) due to ratio $\delta_{cc} \gg 1.0847$
- ▶ Compromise to be able to measure GTR modified asphalt binder
- ▶ Geometry constants must be calibrated
- ▶ Defined as relative measuring systems (same as PP with gap size $> 2 \text{ mm}$ is)
- ▶ Nevertheless, provides data comparable to absolute systems


Calibration of B-CC17SP

Measuring profile used to calibrate and validate the geometry factors:

1	25 Points 0.4 min	<input checked="" type="checkbox"/>	5 Points 1 min	3	25 Points 0.4 min
$\dot{\gamma}$	$\dot{\gamma}$ 1 1/s				
γ				γ 12 %	ω 10 1/s
T	T 76 °C			T 76 °C	

- 1 1/s in rotation was used to be within the newtonian range
- 12 % and 10 rad/s were used according to AAHSTO T315

Verification result after calibration:




➔ No need to perform again at customer site!

➔ Please Note: This calibration does not fit for low viscosity samples!

- Shear viscosity and complex viscosity values are perfectly within $\pm 1\%$ of the viscosity value according to the calibration certificate of the used Cannon N2700000SP
- Limit according to AASHTO T315: $\pm 3\%$

Viscosity Standard: N2700000SP		Lot Number:	
Certification/Issue Date: 06/05/2013		Expiry Date:	
Temperature		Kinematic Viscosity	Dynamic Viscosity
°C	°F	mm ² /s (cSt)	mPa·s (cP)
52.00	125.60		600400
58.00	136.40		395800
64.00	147.20		266700
70.00	158.00		183000
76.00	168.80		127800


Concentric Cylinder Geometry: Reliabe (Real) Data



- ▶ CC compared with data from the accepted plate-plate, 1mm gap, 25mm dia
- ▶ Once thermal equilibrium is achieved, curve shapes are equal
- ▶ No significant difference in visco-elastic behavior (phase angle δ values)
- ▶ Relative CC system with large gap is comparable to absolute PP system (1 mm gap)

Concentric Cylinder Geometry: Proof of Concept Unaged Binder


Influence of GTR mesh size on complex shear modulus G^* (neat asphalt binder)



- ▶ Increasing G^* values with increasing particle size
- ▶ Only slightly difference between 40/80 mesh and 30 mesh
- ▶ Overall stiffness (G^*) is influenced by different rubber particles

Influence of rubber modification on complex shear modulus G^* and phase angle δ




(PAV aged asphalt binder)



- ▶ Increasing stiffness (G^*) with addition of GTR
- ▶ Shift to more solid-like behavior with addition of GTR
- ▶ Dramatic loss of damping ability (flow or liquid-like behavior) at intermediate temperature (delta going towards 0° or solid behavior)

Summary and Conclusions

- ▶ Plate-plate geometry has intrinsic limits:
 - Questionable rheological data due to sample sag
 - Difficult to trim crosslinked rubber particles
 - Appropriate for smaller particle sizes ($\ll 2\text{mm}$)
 - Okay for intermediate temperatures where sample sag is not an issue
- ▶ Concentric cylinder geometries with large gap provide:
 - Elimination of fundamental plate-plate limitations
 - Calibration procedure as for relative measuring systems
 - Easy, straight forward compensation for torsional compliance; thicker bob axis greatly reduces need for compliance compensation
 - **Performance Grading measurements at both high and intermediate temperatures are within reach**
- ▶ Compromises must be made, i.e. use of relative systems
- ▶ Final steps – collaboration partners currently working on repeatability with CC10 bob for PAV aged binder
- ▶ Practical considerations need to be addressed such as longer time to equilibrium, reproducible cup filling, etc...



Some Impressions

